

# Large Diameter Front Loading Priorclaves



- 320 Litre Capacity
- Large diameter chamber giving 20% more loading area than conventional models.
- Quickseal 2 single action door closure system
- Electrically Heated 21kW (three phase 30A)
- Direct steam heated models available
- Steam generated within the chamber (electrically heated models)
- Forced air cooling – standard
- Automatic timed freesteaming – standard
- Conforms to BS2646-1993, BS5500 Cat.3, European PED 97/23/EC, EMC compliance tested and CE Marked.
- Supplied with two full width stainless steel shelves.
- Delayed start and media warming capacity

## Safety Features

The pressure vessel is insurance approved.

Safety features include: protective outer safety cover, special safety handnut mechanism, which cracks the gasket seal allowing any residual pressure to safely escape. A thermal cooling lock, preventing the door from being opened until the contents have cooled to a safe temperature is also fitted. The unit also features a low water level sensor, that prevents boil-dry conditions and a safety blow down valve affords further protection.

## Design

Operation up to 138°C – 2.4 Bar. The **TACTROL** microprocessor control system has staged illumination giving digital readouts of cycle status and the simple to set control panel features adjustable temperature and time that can be viewed simultaneously. Ideally suited for use in a wide range of applications:

*Sterilisation* - media preparation, liquids and diluent, waste, glassware, instruments, apparatus, health care

*Non-sterilisation* – food processing, plastics, accelerated curing, construction.

Extensively used within many industries, including: food, drink, dairy, pharmaceutical, agricultural, education, healthcare, research establishments and industry.

## Installation requirements

The single vent outlet pipe can be directed to a nearby drain. For make safe processing the discharge should be directed to a sealed discharge system; the system should lead by direct connection to a building drain or catchment tank.

*For more specific advice about your particular installation requirements please contact our Technical Department*




## Installation Details

\*where applicable to model chosen

<b>Connection</b>	<b>H320</b>
<b>Electrical*</b>	415V 50Hz 30A per phase TP, E&N
<b>Steam*</b>	Dry Steam reduced to 3.5 Bar. ½" BSP.
<b>Height from floor</b>	1300mm
<b>Vent Outlet</b>	22mm (¾" BSP)
<b>Drain*</b>	15mm (½" BSP)
<b>Water Inlet*</b>	15mm (½" BSP)

Dimensions are approximate only

# Options and Accessories



Priorclaves are available with a full range of control options to suit any application.

- Pre-cycle vacuum stages can be accommodated with vacuum assisted drying and vacuum assisted cooling (patent applied for) with the fitting of a fully automatic vacuum pump
- A selection of baskets and stainless steel discard containers are available for individual models.
- To help with heavy loads basket hoists and loading trolleys are available.

For more specific advice about your specific requirements please contact our Sales Department.

### Load Sensed Process Timing (Option)

Guaranteed sterilisation times from a probe in the anticipated coolest

### Automatic Freesteaming

Includes a venting period during heating, greatly improving temperature distribution and helping air removal. Also available – Pulsed Freesteaming – for improved waste load performance.

### Multi-Program Memory (Option)

Simply press program number required to restore all your settings. Five and ten program options available.

### Accelerated Cooling

Powerful fans blow cold air over the autoclave vessel reducing cooling times to safe load temperatures. Delayed start selectable to reduce media volume loss.

### Media Warming

Forget early morning starts to melt media! The autoclave maintains the media overnight at a 'ready to pour' temperature.

### Vacuum (Options)

Pre-Cycle Vacuum to assist air removal. Vacuum Drying and Cooling options are also available.

### Process Printer (Option)

A full record of temperature and time is recorded at key stages of the cycle. Your details can be printed at the start of each printout, together with the autoclave serial and cycle number. Pressure and load temperature printing options are also available.



## Specifications

(For direct Steam Heated models substitute SH for EH)

Order Ref.	Working Volume (Litres)	Loading Format	Working Dims. wxdxh (mm)	Outside Dims. wxdxh (mm)	Steam Heating Available	Weight (KG)
PS/Q63/EH320	320	Front	630dia x 1020	760x1280x1600	Yes	300



**Priorclave Limited.** 129/131 Nathan Way, Woolwich, London SE28 0AB  
 Tel: +44 (0) 208 316 6620 Fax: +44 (0) 208 855 9734  
 e-mail: [sales@priorclave.co.uk](mailto:sales@priorclave.co.uk) Web: [www.priorclave.co.uk](http://www.priorclave.co.uk)



E. & O.E.