

This authorizes the application of the Certification Mark(s) shown below to the models described in the Product(s) Covered section when made in accordance with the conditions set forth in the Certification Agreement and Listing Report. This authorization also applies to multiple listee model(s) identified on the correlation page of the Listing Report.

This document is the property of Intertek Testing Services and is not transferable. The certification mark(s) may be applied only at the location of the Party Authorized To Apply Mark.

Applicant:	Priorclave Limited	Manufacturer:	Priorclave Limited
Address:	129-131 Nathan Way, West Thamesmead Business Park London SE28 0AB	Address:	129-131 Nathan Way, West Thamesmead Business Park London SE28 0AB
Country:	UNITED KINGDOM	Country:	UNITED KINGDOM
Contact:	Mr Peter Melling / Mr Tony Collins	Contact:	Mr Peter Melling / Mr Tony Collins
Phone:	+44 (0)208 316 6220	Phone:	+44 (0)208 316 6220
FAX:	+44 (0)208 855 0616	FAX:	+44 (0)208 855 0616
Email:	pmelling@priorclave.co.uk / t.collins@priorclave.co.uk	Email:	pmelling@priorclave.co.uk / t.collins@priorclave.co.uk

Party Authorized To Apply Mark: Same as Manufacturer
Report Issuing Office: Leatherhead, UK

Control Number: 4007803

Authorized by: _____


 Ulla-Pia Johansson-Nilsson for
 Thomas J. Patterson, Certification Manager



This document supersedes all previous Authorizations to Mark for the noted Report Number.

This Authorization to Mark is for the exclusive use of Intertek's Client and is provided pursuant to the Certification agreement between Intertek and its Client. Intertek's responsibility and liability are limited to the terms and conditions of the agreement. Intertek assumes no liability to any party, other than to the Client in accordance with the agreement, for any loss, expense or damage occasioned by the use of this Authorization to Mark. Only the Client is authorized to permit copying or distribution of this Authorization to Mark and then only in its entirety. Use of Intertek's Certification mark is restricted to the conditions laid out in the agreement and in this Authorization to Mark. Any further use of the Intertek name for the sale or advertisement of the tested material, product or service must first be approved in writing by Intertek. Initial Factory Assessments and Follow up Services are for the purpose of assuring appropriate usage of the Certification mark in accordance with the agreement, they are not for the purposes of production quality control and do not relieve the Client of their obligations in this respect.

Intertek Testing Services NA Inc.
 545 East Algonquin Road, Arlington Heights, IL 60005
 Telephone 800-345-3851 or 847-439-5667 Fax 312-283-1672

Standard(s):	ANSI/UL 61010-1, 3rd Edition issued May 11, 2012, Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use - Part 1: General Requirements, CAN/CSA-C22.2 No. 61010-1, 3rd Edition issued May 11, 2012, Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use - Part 1: General Requirements
Product:	Range of Autoclaves
Brand Name:	Priorclave

HEG

Models: PNA/aaa/xyz (D)(M)

Where:

aaa represents system type:

QCS is non-vacuum model and QVA is a vacuum model.

x represents the heating type:

E for Electrical Heating, S for Steam Heating

y represents the Vessel orientation:

H = Horizontal, V = Vertical

z represents the Vessel nominal capacity in litres:

100 litres 150 Litres or 200 litres (200 in Horizontal configuration only)

(**D**) represents Double ended system

(**M**) represents optional vacuum drying configuration

** May be followed by Suffixes to indicate configurable options