

Tactrol2 Data Download Sample

Note: Advanced Logging activated.

3590 080213 - Notepad

File Edit Format View Help

Cycle Start	Stage Time	Fault	Chamber	LoadSense	Pressure	FOValue	Cold Junction	LoadSimProbe	HeaterOverTemp	Steam Gen	Pressure	DigitalInputs	D
23:51:57			31/01/13		3590 85	0.000000	15.2	43.5	0.0	3.936	11001111 11111111	0010 00000001 01000000	00000000
PRE-VACUUM	23:52:56			37.6	37.9	-0.109	15.1	45.9	0.0	3.648	11001111 11111111	0010 00111000 01000000	00000000
PRE-VACUUM	23:53:55			60.6	59.6	0.068	15.1	54.7	0.0	3.955	11001111 11111111	0010 00111000 01000000	00000000
PRE-VACUUM	23:54:55			79.7	79.2	0.077	15.1	68.4	0.0	3.582	11001111 11111111	0010 00111000 01000000	00000000
PRE-VACUUM	23:55:54			91.6	91.5	0.000000	15.1	78.5	0.0	3.888	11001111 11111111	0010 00011000 01000000	00000000
HEAT STAGE 2	23:56:38			89.7	98.0	0.013	15.1	81.2	0.0	3.513	11001111 11111111	0010 00011000 01000000	00000000
HEAT STAGE 2	23:56:48			98.0	99.5	0.068	15.1	82.5	0.0	3.582	11001111 11111111	0010 00011000 01000000	00000000
HEAT STAGE 2	23:56:58			99.6	101.1	0.118	15.1	83.9	0.0	3.633	11001011 11111111	0010 00011000 01000000	00000000
HEAT STAGE 2	23:57:07			101.1	102.5	0.174	15.1	85.3	0.0	3.682	11001011 11111111	0010 00011000 01000000	00000000
HEAT STAGE 2	23:57:17			102.5	103.9	0.229	15.1	86.6	0.0	3.730	11001011 11111111	0010 00011000 01000000	00000000
HEAT STAGE 2	23:57:27			103.9	105.3	0.281	15.1	88.0	0.0	3.776	11001011 11111111	0010 00011000 01000000	00000000
HEAT STAGE 2	23:57:37			105.3	106.6	0.340	15.1	89.4	0.0	3.826	11001011 11111111	0010 00011000 01000000	00000000
HEAT STAGE 2	23:57:47			106.6	107.9	0.395	15.1	90.6	0.0	3.872	11001011 11111111	0010 00011000 01000000	00000000
HEAT STAGE 2	23:57:57			107.9	109.2	0.455	15.1	91.9	0.0	3.921	11001011 11111111	0010 00011000 01000000	00000000
HEAT STAGE 2	23:58:07			109.2	110.5	0.514	15.2	93.0	0.0	3.966	11001011 11111111	0010 00011000 01000000	00000000
HEAT STAGE 2	23:58:17			110.5	111.8	0.576	15.2	94.3	0.0	3.780	11001011 11111111	0010 00011000 01000000	00000000
FREE-STEAMING	23:58:18			111.8	112.0	0.586	15.2	94.5	0.0	3.746	11001011 11111111	0010 00111000 01000000	00000000
FREE-STEAMING	23:58:19			112.0	112.0	0.591	15.2	94.5	0.0	3.720	11001011 11111111	0010 00111000 01000000	00000000
FREE-STEAMING	23:58:20			112.0	112.1	0.517	15.2	94.6	0.0	3.703	11001011 11111111	0010 00111000 01000000	00000000
FREE-STEAMING	23:58:21			112.1	112.2	0.470	15.2	94.8	0.0	3.692	11001011 11111111	0010 00111000 01000000	00000000
FREE-STEAMING	23:58:22			112.2	112.0	0.406	15.2	94.9	0.0	3.674	11001011 11111111	0010 00111000 01000000	00000000
FREE-STEAMING	23:58:23			112.0	111.8	0.380	15.2	95.0	0.0	3.664	11001011 11111111	0010 00111000 01000000	00000000
FREE-STEAMING	23:58:24			111.8	111.5	0.356	15.2	95.1	0.0	3.645	11001011 11111111	0010 00111000 01000000	00000000
FREE-STEAMING	23:58:25			111.5	111.1	0.314	15.2	95.2	0.0	3.627	11001011 11111111	0011 00111000 01000000	00000000
FREE-STEAMING	23:58:26			111.1	110.7	0.294	15.2	95.3	0.0	3.619	11001011 11111111	0011 00111000 01000000	00000000
FREE-STEAMING	23:58:27			110.7	110.2	0.277	15.2	95.4	0.0	3.614	11001011 11111111	0011 00111000 01000000	00000000
FREE-STEAMING	23:58:28			110.2	109.7	0.245	15.2	95.5	0.0	3.592	11001011 11111111	0011 00111000 01000000	00001000
FREE-STEAMING	23:58:29			109.7	109.2	0.230	15.2	95.6	0.0	3.585	11001011 11111111	0011 00111000 01000000	00001000
FREE-STEAMING	23:58:30			109.2	108.7	0.203	15.2	95.7	0.0	3.563	11001011 11111111	0011 00111000 01000000	00001000
FREE-STEAMING	23:58:31			108.7	108.2	0.192	15.2	95.8	0.0	3.536	11001011 11111111	0011 00111000 01000000	00001000
FREE-STEAMING	23:58:32			108.2	107.7	0.181	15.2	95.9	0.0	3.522	11001011 11111111	0011 00011000 01000000	00001000
FREE-STEAMING	23:58:33			107.7	107.3	0.197	15.2	95.9	0.0	3.494	11001011 11111111	0011 00011000 01000000	00001000
FREE-STEAMING	23:58:34			107.3	106.9	0.215	15.2	96.0	0.0	3.481	11001011 11111111	0011 00011000 01000000	00001000
FREE-STEAMING	23:58:35			106.9	106.6	0.227	15.2	96.0	0.0	3.472	11001011 11111111	0011 00011000 01000000	00001000
FREE-STEAMING	23:58:36			106.6	106.4	0.247	15.2	96.1	0.0	3.436	11001011 11111111	0011 00011000 01000000	00001000
FREE-STEAMING	23:58:37			106.4	106.3	0.256	15.2	96.1	0.0	3.423	11001011 11111111	0011 00011000 01000000	00001000
FREE-STEAMING	23:58:38			106.3	106.3	0.263	15.2	96.1	0.0	3.398	11001011 11111111	0011 00011000 01000000	00001000
FREE-STEAMING	23:58:39			106.3	107.0	0.276	15.2	96.1	0.0	3.361	11001011 11111111	0011 00011000 01000000	00001000