










Priorclave Illustrated Parts Pricelist




Product Code	Description	Unit of Sales	List Price (£)
<u>Autoclave Consumables</u>			
WTR/FLT/LGE	Replacement Deioniser Cartridge - Large - for Autoclave scale prevention system	Each	£159.26
WTR/FLT/SML	Replacement Deioniser Cartridge - Small for Autoclave scale prevention system	Each	£106.18
ZZZ/BLU/001	Connector Blue Ring M5	Each	£0.30
ZZZ/BLU/002	Connector Blue Sleeved Spade	Each	£0.30
ZZZ/BLU/003	Connector Blue Bullet Female	Each	£0.30
ZZZ/BLU/004	Connector Blue Bullet Male	Each	£0.30
ZZZ/CBT/001	Small Cable Ties	Pk100	£16.35
ZZZ/CBT/002	Large Cable Ties	Pk100	£16.35
ZZZ/CCK/COM	Circular Charts For ABB Commander	PK 50	£24.33
ZZZ/CCK/PEN	ABB Commander Chart Recorder Pen	PK 3	£42.31
ZZZ/DEO/001	Deodorising Capsules Pack - 'Classic' fragrance	Pack	£44.68
ZZZ/DEO/002	Deodorising Capsules Pack - 'Apple' fragrance	Pack	£44.68



Product Code	Description	Unit of Sales	List Price (£)
ZZZ/DEO/003	Deodorising Capsules Pack - 'Citrus' fragrance	Pack	£44.68
ZZZ/FER/BLU	Bootlace Ferrule Connector Blue	Each	£0.05
ZZZ/FER/GRN	Bootlace Ferrule Connector Green	Each	£0.05
ZZZ/FER/GY1	Bootlace Ferrule Connector Grey Single Wire	Each	£0.05
ZZZ/FER/GY2	Bootlace Ferrule Connector Grey Double Wire	Each	£0.08
ZZZ/FER/ORA	Bootlace Ferrule Connector Orange	Each	£0.05
ZZZ/FER/WHT	Bootlace Ferrule Connector White	Each	£0.05
ZZZ/RED/001	Connector Red Ring M4	Each	£0.30
ZZZ/RED/002	Connector Red Ring M6	Each	£0.30
ZZZ/RED/003	Connector Red Unsleeved Small Spade	Each	£0.30
ZZZ/RED/004	Connector Red Sleeved Standard Spade	Each	£0.30
ZZZ/RED/005	Connector Red Bullet Female	Each	£0.30
ZZZ/RED/006	Connector Red Bullet Male	Each	£0.30






Product Code	Description	Unit of Sales	List Price (£)
 ZZZ/RIB/001	Printer Ribbon Type 1 (Each)	Each	£13.59
 ZZZ/RIB/002	Printer Ribbon Type 2 (Each)	Each	£13.59
 ZZZ/ROL/005	Printer Paper Roll (Pack of 5)	Pk 5	£13.59
ZZZ/SCK/FAN	P100/41/7414X ABB Fan Fold Chart Paper for CR100 chart recorder	PK2	£28.44
ZZZ/SCK/PNG	Kent Strip Chart Rec. CR100 Green Pen	PK 3	£38.19
ZZZ/SCK/PNM	PR100/021/0 ABB 3-Colour Pen Capsule for SR100 chart recorder	Each	£67.14
ZZZ/SCK/PNR	Kent Strip Chart Rec. CR100 Red Pen	PK 3	£38.19
ZZZ/SPD/001	Small Sticky Pads	Pk 100	£12.80
ZZZ/SPD/002	Large Sticky Pads	Pk 100	£12.80
ZZZ/SYP/LGE	Long Siphon Pump for Top Loading Autoclaves	Each	£27.40
ZZZ/SYP/SML	Short Siphon Pump (for front loading autoclaves)	Each	£15.06
ZZZ/TCG/001	Tub of copper grease	Each	£18.73
ZZZ/TKC/999	Complete Tool Kit	Kit	£330.74

Product Code	Description	Unit of Sales	List Price (£)
ZZZ/TSG/001	Tube Silicone Grease	Each	£21.53
ZZZ/TSS/001	Silicone Sealant - clear	Each	£16.43
ZZZ/TSS/002	Silicone Sealant - Black	Each	£16.43
ZZZ/TSS/003	Silicone Sealant - White	Each	£16.43
ZZZ/TSS/004	Silicone Sealant - Grey	Each	£16.43
ZZZ/YEL/001	Connector Yellow Ring M6	Each	£0.30
ZZZ/YEL/002	Connector Yellow Sleeved Spade	Each	£0.30
<u>Common Spares</u>			
ADT/ACS/15F	15mm x 1/2in F.I.C Adaptor	Each	£8.53
ADT/ACS/22F	22mm x 3/4in F.I.C Adaptor	Each	£8.53
	AFL/ACS/LGE Replacement Air Filter Large (Intake)	Each	£129.29
	AFL/ACS/SML Replacement Air Filter Small (Intake)	Each	£86.18
BFI/ACS/015	Boiler Fitting Complete 15mm (1/2in Tank Connector)	Each	£10.78
BFI/ACS/022	Boiler Fitting Complete 22mm (3/4in Tank Connector)	Each	£14.38

Product Code	Description	Unit of Sales	List Price (£)
FLO/WAT/SWT	Float Switch Header / Waterfill Tank	Each	£29.73
FUS/ACS/006	Legrand Fuse 6A	Each	£0.77
FUS/ACS/016	Legrand Fuse 16A	Each	£0.77
GAS/STG/001	Gasket for Steam Generator 100L	Each	£39.35
HBA/WHE/REP	Replacement Handle for Original Bolt	Each	£81.64
 HCO/ACS/150	150C Heater Cut-out (Thermal Fuse)	Each	£3.57
 HCO/STG/180	Thermal Cut-out 180C For Steam Generators	Each	£3.88
INL/ACS/FIL	Inline Electrical Filter (Power Line Filter)	Each	£32.35
 KEY/SLK/001	Setting Lock Keyswitch & Key Tactrol	Each	£37.31
KEY/SLK/002	Programmer Setting Lock Key EB0001	Each	£22.33
KEY/TLK/004	Thermal Lock Key EAO 311	Each	£5.13
 LWP/ACS/00C	Type 'C' Low Water Electrode	Each	£77.28
LWP/DPB/000	Dual low water probe for boilers. <i>Check electrode details. DSS - Dual, Straight, Straight.</i>	Each	£81.41
LWP/RSC/000	Low Water Probe (no Tag) for Front Loading Priorclaves inc. C40 & RSC	Each	£77.29





Product Code	Description	Unit of Sales	List Price (£)
LWP/ULT/SEN	Ultrasonic Low Water Sensor	Each	£224.88
MCB/VAC/003	Three Phase MCB for Vac Pumps & Motors	Each	£79.56
MSA/LLK/MID	Lid Locking Pin Microswitch Assembly - Compacts	Each	£8.62
 PRB/ACS/003	Replacement Control Probe Duplex <i>Stainless Steel Fitting PRB/ACS/FSS normally required</i>	Each	£73.44
 PRB/ACS/004	Replacement Wandering Probe (Tactrol Type) <i>Requires conversion kit if flat pin plug fitted</i>	Each	£74.95
 PRB/ACS/005	Load Simulator Thermocouple	Each	£33.65
PRB/ACS/FSS	3mm x 1/8in BSP fitting for Duplex Probe S/S <i>Usually required when re-fitting Duplex Probes</i>	Each	£5.14
PRB/PLU/SOC	Thermocouple Plug & Socket for wandering Probe	Each	£4.03
PSS/ADT/001	Pressure Switch Fitting	Each	£7.62
 PSS/SWI/001	Pressure Switch	Each	£75.00
 PST/ACS/001	Pressure Transducer	Each	£99.94
RLY/PLG/008	8 Pin Plug	Each	£15.67
 RLY/SSR/COV	Solid State Relay Cover	Each	£7.77

	Product Code	Description	Unit of Sales	List Price (£)
	SBG/ACS/015	Silicone Bung 15mm	Each	£7.19
	SBG/ACS/022	Silicone Bung 22mm	Each	£10.06
	SPG/MID/SOL	Locking Solenoid Return Spring	Each	£4.87
	STR/ACS/015	15mm (1/2in BSP) Strainer	Each	£20.38
	STR/ACS/028	28mm (1" BSP) Strainer	Each	£31.54
	STR/ACS/WHS	Conical Fitting Strainer Small	Each	£20.72
	TAC/OPT/001	Tactrol Output Board	Each	£147.97
	TAC/PCB/001	Tactrol Main P.C.B	Each	£1,141.59
	TAC/PCB/UPG	Tactrol Software Upgrade	Each	£74.41
	TAC/POW/001	Switch Mode Power Supply	Each	£165.38
	TAC/PRN/001	Tactrol Printer Replacement	Each	£594.14
	TAC/RIB/DIS	Display to Main Board Ribbon Cable <i>Length Serial No. Dependant</i>	Each	£14.81
	TC1/CON/001	Main Control Panel	Each	£339.45
	TC1/PGR/010	□10 Program Memory Panel	Each	£107.39

Product Code	Description	Unit of Sales	List Price (£)
TC1/PRG/005	5 Program Memory (Tactrol 1)	Each	£164.91
TC1/PRN/CON	<input type="checkbox"/> Printer Conversion Panel	Each	£107.39
 TC2/CON/001	Tactrol 2 Main Control Panel	Each	£339.45
 TC2/CON/EXG	Exchange Control Panel (Tactrol 2)	Each	£339.45
TC2/PRG/005	Tactrol 2 Program Select Panel	Each	£164.91
 VCE/ACS/001	Condensate Trap 15mm (1/2in)	Each	£128.51
 VPC/ACS/001	Air Purge Valve Complete	Each	£128.51
 VPE/ACS/002	Air Purge Valve Element (J-Fill) <i>Replaces 'G' fill elements previously used</i>	Each	£97.73
 VSC/ACS/015	Safety Valve Complete 1/2" (15mm)	Each	£180.77
 VVC/ACS/022	Vent Valve Complete 3/4" (19mm)	Each	£110.09
 VVC/RSC/028	28mm Solenoid (Vent) Valve for RSC Priorclaves	Each	£176.67
VVF/ACS/015	15mm Vent Valve Fitting	Each	£3.42
VVF/ACS/022	22mm Vent Valve Fitting	Each	£3.42
VVF/ACS/028	28mm Vent Valve Fitting	Each	£9.87

Product Code	Description	Unit of Sales	List Price (£)
WAS/FIB/012	12mm Fibre Washer for heaters	Each	£0.61
WAS/FIB/015	Fibre Washer- 15mm Tank Connector Washer	Each	£0.77
WAS/FIB/022	Fibre Washer - 22mm Tank Connector	Each	£1.15
WAS/SET/001	Complete set of washers, with PTFE tape & jointing compound. <i>Contents dependant on Serial No.</i>	Each	£38.68
WTR/FLT/CTG	Deionising Water Filter	Each	£100.62




Exchange Components

	TAC/PCB/EXG	Exchange Tactrol Processor Board <i>Exchange involving component replacement</i>	Each	£703.11
	TAC/PCB/RES	Exchange Tactrol Processor Board (Reset) <i>Board reset only</i>	Each	£74.41
	TAC/PRN/EXG	Tactrol 2 Printer Exchange <i>Note; Exchange rarely applicable due to improved lifespan.</i>	Each	£138.38
	TC1/CON/EXG	Tactrol Control Panel Exchange	Each	£202.35
	VAC/ACS/EXG	Standard vacuum pump, reconditioned, 6 month warranty, old unit must be in mechanically good condition. <i>Subject to availability, incl. new start capacitor. Old unit must be in mechanically good condition.</i>	Each	£586.74

Spares - Compact Priorclaves




	CAS/MID/000	Set of 4 castors for Midas V55, V56 and Compact 60	Each	£20.54
	COV/DOR/C36	Midas C32/36 Door Cover	Each	£246.63




Product Code	Description	Unit of Sales	List Price (£)
FAN/MID/C56	Cooling Fan Circular for Midas 52/55/56	Each	£81.50
FAN/TAN/C60	C60 Tangential Cooling Fan	Each	£61.33
 GAS/MID/300	300mm Midas Gasket for Compact 32/36 & 55/56 models	Each	£54.61
 GAS/MID/350	350mm Gasket for Compact 40/60 models	Each	£54.61
HND/BLK/C60	Locking Handle For Compact Priorclaves	Each	£13.40
 HTR/MID/C40	Compact 32/36/40 Heater	Each	£93.69
HTR/MID/C60	Compact V55/56/60 Heater	Each	£93.69
LLA/MID/001	Midas Lid Locking Arm	Each	£73.83
LLN/MID/000	Stainless steel lid liner for Midas autoclaves	Each	£44.64
MPS/MID/001	Compact Mains Power Switch	Each	£28.44
MSA/MID/SOL	Solenoid Microswitch Assembly (Midas) (32/36/55/56)	Each	£8.62
PSG/GAU/063	63mm Pressure Gauge	Each	£39.82
PSS/CON/001	Pressure Switch Conversion Kit <i>Includes Pressure Switch Fitting</i>	Kit	£39.30

Product Code	Description	Unit of Sales	List Price (£)
 RLY/MID/HTR	Heater Switching Relay	Each	£23.86
 RLY/SSR/240	240V Solid State Relay	Each	£58.44
SHM/MID/HNG	Hinge Shim for Compact Priorclaves <i>All Thicknesses & Materials</i>	Each	£0.86
 SOL/MID/001	Compact Lid Locking Solenoid	Each	£41.37
SPG/LID/C60	Lid Assist Spring C60 Priorclaves	Each	£51.03
TER/FUS/COM	Fused Terminal Block Top Loading Compacts	Each	£11.44
THW/MID/001	Load Simulator Thermo Well - Benchtop Compact Priorclaves	Each	£56.04
VAC/DIA/SGL	C40 Single Head Diaphragm Vacuum Pump	Each	£796.99
VAC/MID/RPR	Vacuum Pump Repair Kit Single Head Gast Pumps <i>Repair kit Gast K294A, diaphragm, reed valve and filter.</i>	Kit	£22.42
VVC/ACS/015	Vent Valve Complete 1/2" (16mm dia)	Each	£106.17
<u>Spares - Large Priorclaves</u>			
JGA/AFL/EXH	Inlet Steam Filter JG type Autoclaves	Each	£843.58
JGA/BGC/010	3/8" Brass Gauge Cock	Each	£21.59
JGA/FHS/DOM	Domnic Hunter Filter Housing Seal	PK 5	£97.32





Product Code	Description	Unit of Sales	List Price (£)
JGA/GAS/001	Door Gasket (John Guest A/C) <i>Requires Making Onsite</i>	Each	£138.89
JGA/HND/001	Handwheels (John Guest A/C)	Each	£198.43
JGA/PIL/VLV	Pilot Valve For SMC Type Pneumatic Valve	Each	£62.17
JGA/SCB/KIT	FT14 Steam Condensate Ball Valve Maintenance Kit	Each	£125.96
JGA/SMC/VSI	SMC Pneumatic 3/4 Steam Inlet Valve	Each	£167.96
JGA/SPI/COM	DP 17 LC Spirax Valve Complete	Each	£546.57
JGA/SPI/DIA	Main Diaphragm	Each	£18.37
JGA/SPI/MKT	Spirax Maintenance Kit	Each	£174.96
JGA/SPI/MPR	Main Push Rod And Plate	Each	£112.78
JGA/SPI/MVA	Main Valve Assembly <i>Contains Main Valve Head And Seat</i>	Each	£43.01
JGA/SPI/MVS	Main Valve Return Spring	Each	£4.49
JGA/SPI/PIL	Pilot Diaphragm	Each	£8.25
<u>Spares - Status Priorclaves</u>			
FAN/STA/001	Cooling Fan for Status Autoclaves	Each	£58.99






Product Code	Description	Unit of Sales	List Price (£)
GAS/STA/LID	Status Lid Gasket	Each	£78.91
MCB/STA/005	Resettable Fuse 5A	Each	£29.84
MCB/STA/016	Resettable Fuse 16A	Each	£29.84
RHB/STA/001	Status Red Handbolt	Each	£252.89
SOL/SLI/STA	Status Locking Plate	Each	£57.20
<u>Recorder & Printer Spares</u>			
TAC/RIB/PRN	Ribbon Cable for Tactrol Printer <i>Length Serial No. Dependant</i>	Each	£9.66
<u>Spares - Standard, Quickseal, RSC & EX Series Priorclaves</u>			
AFL/ACS/EXH	Exhaust Filter (Containment type)	Each	£453.81
AFL/ACS/EXI	Internal S/S Exhaust Filter, 0.2 micron PTFE	Each	£407.69
BLT/STG/000	Steam generator end plate bolt. High tensile M10 x 40. Set of 16.	set	£27.29
CAN/ACS/001	Replacement Vertical Canopy for Top Loading Priorclaves	Each	£558.05
	CON/240/09A 240V Contactor - 9 Amp <i>Tactrol only:Send with output board. FOC if board returned</i>	Each	£41.36
	CON/240/23A 240V Contactor - 23 Amp	Each	£67.82
	CON/240/45A 240V Contactor - 45 Amp	Each	£93.44





Product Code	Description	Unit of Sales	List Price (£)
DTM/240/002	2 Contact Variable Delay Timer For Delaying Process Timer	Each	£71.86
FAN/ACS/001	Replacement Cooling Fan for Top & Front Loaders	Each	£82.16
FAN/ACS/080	80x80x35mm Cooling Fan	Each	£38.00
FAN/ACS/120	120x120x38mm Cooling Fan	Each	£42.61
FAN/ACS/PLG	Plug for Cooling Fan (3 Pin)	Each	£8.23
FAN/QCS/000	Cooling Fan QCS Series Priorclaves	Each	£103.57
FIL/SUP/001	Capacitor Suppressor 0.1µf	Each	£8.19
FIL/VDR/002	V275 LA4 Varistor (VDR)	Each	£2.12
 FLS/ACS/316	S/S 316 Float Switch	Each	£129.37
FLS/ACS/TNK	Float Switch for Vaccum Pump Tank	Each	£31.85
 GAS/ACS/500	500mm Dia Priorclave Gasket	Each	£78.91
 GAS/ACS/660	660mm Dia Priorclave Gasket	Each	£96.58
GAS/ACS/Q63	Priorclave Gasket for Q63 models	Each	£87.50
GAS/ACS/WEI	Red Weir Seal Gasket	Each	£9.32



Product Code	Description	Unit of Sales	List Price (£)
GAS/MID/WEI	Red Weir Seal Gasket Midas 32	Each	£7.58
 GAS/RSC/230	Priorclave Gasket RSC 160/230/350L Models	Each	£109.39
 GAS/RSC/450	Priorclave Gasket RSC 450/700 Models	Each	£124.24
 GAS/RSC/SPC	Gasket for 1500L Priorclave	Each	£137.61
HBA/RSC/000	Number Discontinued Now HBA/RSC/003	Each	£0.00
HBA/RSC/001	RSC Red Safety Bolt Handwheel Assembly	Each	£212.94
HBA/RSC/002	RSC Plain Bolt Handwheel Assembly	Each	£206.61
HBA/RSC/003	RSC Thermal Lock Closure Bolt Handwheel Assembly	Each	£206.61
HBA/RSC/010	Number Discontinued Now HBA/RSC/001	Each	£0.00
HBA/RSC/020	Number Discontinued Now HBA/RSC/002	Each	£0.00
HBA/RSC/AST	RSC Handwheel Assist	Each	£32.43
HBB/ACS/001	Black Handnut Assembly with insert	Each	£159.25
HBR/ACS/001	Red Handnut Assembly with S/S Insert	Each	£252.89




Product Code	Description	Unit of Sales	List Price (£)
HND/MID/000	Replacement Handle for Quickseal Safety Gate	Each	£20.54
HND/QCS/YEL	Locking Handle For QCS	Each	£28.90
HNG/ACS/001	Replacement door hinge set inc. hinge pins, sliders, thrust washer	Each	£190.48
HPD/CBL/001	Power door lifting cable. Right hand side	Each	£37.31
HPD/CBL/002	Power door lifting cable. Left hand side.	Each	£35.88
 HTR/ACS/HOR	Heating Element 3.5kW. Front Loading Priorclaves	Each	£93.69
 HTR/ACS/VER	Heating Element 3.5kW. Top Loading Priorclaves	Each	£93.69
HTR/ETC/000	Heating element for environmental chamber	Each	£26.44
HYD/COI/BLK	Black Hydraulic Coil & Plug 220V	Each	£30.92
HYD/COI/RED	Red Hydraulic Coil & Plug 220V	Each	£35.95
HYD/CON/PLC	Hydraulic Door PLC	Each	£541.63
HYD/PUL/BKT	Hydraulic Pulley Wheel Bracket	Each	£58.17
HYD/PUL/WHE	Hydraulic Pulley Wheel	Each	£56.78




Product Code	Description	Unit of Sales	List Price (£)
HYD/SOL/001	Hydraulic Solenoid Valve Without Manual Override	Each	£183.62
HYD/SOL/002	Hydraulic Solenoid Valve With Manual Override	Each	£201.22
KEY/SLK/003	EAO 2300 Setting Lock Key	Each	£22.33
KEY/STP/455	455 Key For Emergency Stop Button	PK 2	£12.74
KEY/TLK/003	Thermal Lock Key ED9001	Each	£5.13
LID/DOR/RPR	Lid /Door repair including fitting liner	Each	£1,556.39
LLN/ACS/000	Lid Liner for Standard Priorclaves	Each	£111.62
 MCB/ACS/006	6 Amp Main Circuit Breaker	Each	£26.53
 MCB/ACS/016	16Amp Main Circuit Breaker (triple pole)	Each	£79.56
 MCB/ACS/032	32Amp Main Circuit Breaker (triple pole)	Each	£119.57
 MCB/ACS/045	45Amp Main Circuit Breaker (triple pole)	Each	£150.93
MCB/VAC/001	Single Phase MCB for Vac Pumps and Motors	Each	£26.52
MOD/QCS/SPG	Spring Kit To Replace QCS Gas Spring	Kit	£74.84






Product Code	Description	Unit of Sales	List Price (£)
MSA/CAN/HOR	Canopy Microswitch Assembly - horizontal	Each	£8.62
MSA/CAN/VER	Canopy Microswitch Assembly	Each	£8.62
	MSA/DCL/QCS Door Closed Microswitch QCS Type Priorclaves	Each	£47.89
MSA/LID/QCS	QCS Lid/Door Microswitch	Each	£8.62
MSA/QCK/SBT	Midas Quickseal Safety Bolt Microswitch	Each	£23.61
	MSA/RSC/ROL Roller Microswitch Assy. Rectangular Series	Each	£28.74
	MSA/RSC/SBT Safety Bolt M/Switch Assy. Rectangular Series	Each	£43.12
MSA/SBT/ACS	Safety Bolt Lever Microswitch Assembly	Each	£11.81
MSA/WAT/ACS	Waterproof Type Safety Bolt M/Switch Assembly	Each	£42.78
MSA/WAT/PRN	Waterproof Type Safety Bolt M/Switch Assembly (Printer)	Each	£53.92
NUT/STG/000	Steam generator nut set. M10 full nut. Set of 16	set	£7.68
	PSG/ACS/100 Vacuum / Pressure Gauge 100mm -1 to 3 Bar	Each	£126.14
	RLY/008/12V 12V 8 Pin Relay	Each	£24.16



Product Code	Description	Unit of Sales	List Price (£)
 RLY/SSR/DC1	Solid State Relay DC/AC	Each	£76.61
 SHM/QCS/CLA	Locking Clamp Shim QCS Priorclaves <i>All Thicknesses & Materials</i>	Each	£1.29
SHM/QCS/HNG	Hinge Block Shim QCS Priorclaves <i>All Thicknesses & Materials</i>	Each	£1.29
SHM/RSC/001	Hinge Block Shim for RSC Priorclaves <i>All Thicknesses & Materials</i>	Each	£1.70
SLI/QCS/HND	QCS slider bar to link handle to actuator	Each	£69.13
SOL/ACS/001	Lid Locking Solenoid	Each	£106.55
 SOL/QCS/001	Locking Solenoid QCS Series Priorclaves	Each	£32.68
 SOL/RSC/001	Door Locking Solenoid RSC Type	Each	£136.56
SOL/SLI/HOR	Horizontal Locking Solenoid Slider Plate	Each	£26.43
SOL/SLI/LFT	Vertical Locking Slider Left (From Front)	Each	£52.86
SOL/SLI/RHT	Vertical Locking Slider Right Hand (From Front)	Each	£52.86
SPG/QCS/GAS	Gas Spring Quickseal & QCS Models	Each	£101.47
STP/ACS/001	Emergency Stop Button (W/Out Contact Blocks)	Each	£54.96

Product Code	Description	Unit of Sales	List Price (£)
STP/ACS/002	Emergency Stop Button Contact Block	Set of	£23.56
	STR/ACS/022 22mm (3/4in BSP) Strainer	Each	£25.53
STR/ACS/WHL	Conical Fitting Strainer Large	Each	£31.42
TAC/PCB/CHG	Charger, to keep spare mainboard battery charged when not in use <i>Non stock, purchase and fit connector. Check mains plug type for customer.</i>	Each	£87.25
TAC/POW/002	Power Supply For QCS Locking Solenoid	Each	£154.57
TER/ACS/001	Neutral Terminal Block (Rail Mounting)	Each	£3.22
	THW/ACS/001 Load Simulator Thermo Well old type	Each	£72.05
THW/ACS/002	Standart Load Simulator Thermowell	Each	£59.44
TMP/ACS/HSG	Jumo 'Itron 8' Temperature Controler Housing	Each	£48.23
VAC/ACS/RD1	Sequential reed valve for Zeta VP10/4 vacuum pump	Each	£8.07
VAC/ACS/RPR	Pair of Parts kits to repair vacuum pump. <i>Pair of repair kits Gast K294A, diaphragm, reed valves (2), and filters(2).</i>	Pair	£43.72
VAC/DIA/DBL	Twin Head Diaphragm Vacuum Pump	Each	£821.61
VAC/LRP/LRG	Liquid Ring Vacuum Pump 450 -700L 3 Phase <i>1.5kw</i>	Each	£1,390.93

Product Code	Description	Unit of Sales	List Price (£)
VAC/LRP/MED	Liquid Ring Vacuum Pump 200 - 350L 3 Phase .75kw	Each	£1,237.18
VAC/LRP/SM1	Liquid Ring Vacuum Pump 100 -150L 1 Phase 0.35kw	Each	£1,077.28
VAC/LRP/SM3	Liquid Ring Vacuum Pump 100 -150L 3 Phase 0.4kw	Each	£1,022.95
VAC/TUB/RPR	Replacement Tubes, Fittings and Olives- Large Gast Vacuum Pumps	Kit	£22.04
VCE/ACS/002	Condensate Trap Element (E-Fill)	Each	£97.73
VFL/ACS/TOR	Water Fill Float Valve (Torbeck)	Each	£25.32
 VNR/ACS/015	15mm (1/2in) Non Return Valve	Each	£24.54
 VNR/ACS/022	22mm (3/4in) Non Return Valve	Each	£27.99
VNR/ACS/028	28mm (1in) Non-Return Valve	Each	£30.54
 VSC/ACS/022	Safety Valve Complete 3/4" (22mm)	Each	£232.65
VSC/STG/015	Steam Generator Safety Velve 4.5Bar 15mm	Each	£243.44
VSI/ACS/015	Steam Inlet Valve	Each	£116.62
WTR/DSC/ACS	Product used to descale autoclaves	Each	£13.36

Product Code	Description	Unit of Sales	List Price (£)
WTR/DSC/CON	Price to descale autoclave and chemicals during STD visit	Each	£189.70
Spares Pre- Tactrol Models			
ANC/ACS/001	Amber Neon "COOL"	Each	£13.41
ANH/ACS/001	Amber Neon "HEAT"	Each	£13.41
	BUL/ACS/010 Pack 10 Bulbs	Pack 1	£28.11
	CON/24V/09A 24V Contactor - 9 Amp	Each	£41.36
	DTM/240/004 4 Contact Variable Delay Timer (has contacts A1 & A2)	Each	£54.99
	DTM/24V/003 3s 24V Delay Timer	Each	£46.53
	FUS/ACS/010 Pack 10 Assorted Fuses 6A x 4, 16A x 6.	Pack 1	£22.40
	HNG/PIN/V55 Hinge Pin Set V55 Control Panel Perspex Cover	Pair	£25.22
	KEY/TLK/001 Thermal Lock Keyswitch and Key Tactrol	Each	£56.20
	KEY/TLK/002 Thermal Lock Keyswitch & Key (TOK2)	Each	£37.31
	LEN/ACS/CLE Clear Lens for Switch/Indicators <i>State Text Required when ordering</i>	Each	£3.05
	LEN/ACS/GRE Green Lens for Switch/Indicators <i>State Text Required when ordering</i>	Each	£3.05

Product Code	Description	Unit of Sales	List Price (£)
LEN/ACS/ORA	Orange Lens for Switch/Indicators <i>State Text Required when ordering</i>	Each	£3.05
LEN/ACS/RED	Red Lens for Switch/Indicators <i>State Text Required when ordering</i>	Each	£3.05
 LEV/ACS/OMR	IMO type Floatless Level Switch	Each	£88.11
LEV/ACS/SYR	Syrelec Floatless Level Switch	Each	£88.11
 LIU/ACS/001	Lamp Indicator Unit	Each	£17.51
 PBS/ACS/001	Push Button Switch	Each	£45.39
PRB/ACS/002	Replacement Control Probe Single	Each	£70.31
PRB/ACS/FIT	Thermocouple 3/8" fitting <i>For Thermoregulator Replacement</i>	Each	£7.90
PRB/CHR/001	Replacement Wandering Probe (S/S Type) <i>Chart Recorders Only</i>	Each	£78.42
 RLY/240/011	240V 11 Pin Relay	Each	£24.16
 RLY/BAS/811	11 Pin Relay Base Screw Type (8-11 pin converter)	Each	£22.28
RLY/PLG/011	11 pin Plug	Each	£15.67
 STF/ACS/001	Freesteam Thermostat	Each	£124.62
 STT/ACS/001	Thermal Lock Thermostat	Each	£157.44

Product Code	Description	Unit of Sales	List Price (£)
SWV/ACS/COV	Vent Switch Cover	Each	£4.24
TER/FUS/001	Fuse Terminal Block	Each	£11.44
 TIM/ACS/ANL	Analog Timer	Each	£114.87
TIM/ACS/DIG	Digital Timer	Each	£250.65
 TMP/ACS/001	Temperature Controller	Each	£430.90
TMP/ACS/002	Temperature Controller Conversion Kit (without brazed nut)	Kit	£30.55
TMP/LSC/OMR	Omron E5CN - R2MTC Loadsense Controller	Each	£255.22
<u>Service Packs</u>			
SERV999	Service Introduction Pack	Pack	£0.00